

Cylindrical Toothbrush “TAMPOPO”



Brushing teeth is, as we all know, very important for health and beauty.

Shapes of teeth with grooves are uneven so it is difficult to be brushed.

Tongue care became practiced by the people who are conscious of breath odor.

It is remarkable, not only its aspect but also various functions.

We studied with a doctor about a head of toothbrush in consideration of shapes of teeth. The head shape we found is simple cylindrical.

It would fit front teeth or back molar neatly also it titillates gums while brushing.

With a conventional toothbrush, some people complain that the backside of front teeth would be brushed perfunctory. No wonder, it is awkward for such area by structurally.

Using the Cylindrical Toothbrush, you can feel their squeaky-clean surface after brushing. The Cylindrical Head fits that area neatly so that you can brush there easily in back-and-forth motion or in up-and-down motion approaching from backside. The ultrafine bristles are gentle enough for tongue cleaning or gum massaging.

We would like you to consider trying this multi-purpose Cylindrical Toothbrush.

<http://www.cylindrical-toothbrush.info/>
info@toothbrush.co.jp

STB Higuchi

Cylindrical Toothbrush “TAMPOPO”

Specification

Bristle

Material: Nylon
Bristle diameter: approx. 0.09mm
Number of bristles: approx. 21600 1200 x18 layers
Head size in diameter: approx. 18mm

-Reference- Conventional Toothbrush
Thread thickness: 0.2mm in average
Number of bristles: approx.700 to 1000
Head size in height: 15 to 17mm

Handle

Material: Polypropylene
Length: approx. 175mm
Width in diameter: 12mm max. 5mm min.
Color: Pink, Blue, Orange



Features

Ultrafine HD bristles remove plaque powerfully.

It titillates gums while brushing teeth.

It works as a high performance tongue cleaner as well.

Cylindrical head fits the backside of front teeth neatly so that you can brush there easily.

Both upper and lower teeth can be brushed at a time, when you familiarized with it.

All direction bristle head is best suited for Elder Care.

For kids or disabled people who is awkward using a conventional toothbrush.

Contact

STB Higuchi co.,ltd. 81 6 6730 1712
info@toothbrush.co.jp
Mr. Koichi Ishii